

## Medical Equipment

Technical Data Sheet

# Spackman Cannula

The Spackman Cannula is unique in having both a **moveable cervical cone** and an **atraumatic elastomeric tip**.

The **moveable cervical cone** allows the Cannula to be 'made to fit' the depth of the uterus which ensures an effective seal and reduces the risk of dye leakage. The atraumatic elastomeric tip is designed to reduce the risk of perforation to the uterus.


The Spackman Cannula can be purchased separately or as part of the DiOne procedure pack.

## DiOne

The DiOne is a completely disposable procedure pack, containing the equipment required to perform a Laparoscopy and Dye procedure in an easy to access disposable tray.

The DiOne contains:

- a Spackman Cannula, with an atraumatic elastomeric tip, with a **moveable cervical cone** ensuring an effective seal.
- a pair of disposable, steel Teal's Uterine Tenaculum, which aid to attach to the cervix and work in conjunction with the Spackman Cannula to create a tight seal.
- a Clearprobe Disposable Sound Dilator (DSD). The DSD is a flexible, moulded plastic, double-ended sound dilator for assessing the depth of the uterus and to allow for gentle dilation of the cervix.



The Spackman Cannula's moveable cervical cone.

### Ordering Information

#### DiOne kit

Part Number	DI-1
Quantity	5 units per case
Description	Lap and Dye kit
Delivery	Usually ex-stock

#### DSD

Part Number	DSD-1000
Quantity	50 units per case
Description	Disposable Sound Dilator
Delivery	Usually ex-stock

#### Spackman Cannula

Part Number	SPK-01
Quantity	10 units per case
Description	Spackman Cannula
Delivery	Usually ex-stock



The DiOne disposable Lap and Dye kit.

[www.clinicalinnovationseurope.com](http://www.clinicalinnovationseurope.com)

Distributed by