


The Kiwi™ Omnicup Vacuum Delivery System

The Kiwi™ OmniCup is a complete vacuum delivery system.

The cup has a hyper-low Mushroom shaped profile, which is suitable for all OA, OT AND OP fetal malpositions and allows for placement over the flexionpoint more frequently than any other cup.

The PalmPump™ provides a safe and effective source of vacuum giving the operator full control. the unique traction device allows a considerable amount of flexibility in manipulation of the fetal head.

The Kiwi system is easy-to-use, effective, quiet and has proven to be less psychologically threatening to the mother-to-be.



The flexion point is situated on the sagittal suture line 3cm forward of the posterior fontanelle.


Ordering Information

Part Number	VAC-6000M
Quantity	5 units per case
Delivery	Usually ex-stock
Description	Kiwi™ Omnicup Vacuum Delivery System

Educational Equipment

	Part No.
Lucy instrumental birthing simulator	OT 003
Pelvic model with articulated head	PFM-6000
Fetal head without scalp	PFM-6000F
PFM-6000 plus 1 each VAC-6000M, S and MT	PFM-6000C
Handbook of Vacuum Extraction	HBVD
CD ROM Pack (Reducing the Risk & Kiwi™ C-section)	CDR

Kiwi is the trademark of Clinical Innovations, Inc



Distributed