


# The Kiwi™ ProCup Vacuum Delivery System

**The Kiwi™ ProCup is a complete vacuum delivery system designed for comprehensive control without the need of an assistant or any other equipment, such as a vacuum pump.**

The cup has a bell shaped which is suitable, in general, for all positions where the Flexion point is easily available in OA presentations

The PalmPump™ provides a safe and effective source of vacuum right in the palm of your hand. The operator has full control and the unique traction device allows a considerable amount of flexibility in manipulation.

The Kiwi system is easy-to-use, effective, quiet and has proven to be less psychologically threatening to the mother-to-be.



The flexion point is situated on the sagittal suture line 3cm forward of the posterior fontanelle.

## Ordering Information


Part Number	VAC-6000S
Quantity	5 units per case
Delivery	Usually ex-stock

Description	Kiwi™ Procup Vacuum Delivery System
-------------	-------------------------------------

## Educational Equipment

	Part No.
Lucy instrumental birthing simulator	OT 003
Pelvic model with articulated head	PFM-6000
Fetal head without scalp	PFM-6000F
PFM-6000 plus 1 each VAC-6000M, S and MT	PFM-6000C
Handbook of Vacuum Extraction	HBVD
CD ROM Pack (Reducing the Risk & Kiwi™ C-section)	CDR

Kiwi is the trademark of Clinical Innovations, Inc



[www.clinicalinnovationseurope.com](http://www.clinicalinnovationseurope.com)

Distributed